

Overview

The purpose of the program is to collect statistics about the localization of mRNAs in smFISH images of neurons. The program is broken into two parts. The first part, Part 1, is a semi-automatic tool for segmenting the cellular compartments of the neuron, the somas, nuclei, and dendrites. Each segmentation is stored as a binary image which we refer to as a print. Additionally, Part 1 generates skeletons for each dendrite, which are also stored as binary images. After running Part 1, FISHQuant is used to identify the coordinates of mRNAs (x, y, and z) in the smFISH images. Then the second part of the program, Part 2, uses both the spot coordinate (x and y) obtained from FISHQuant along with the skeletons and prints obtained from Part 1 to collect statistics about the localization of mRNA. In Part 2, we investigate the following: 1) The density of mRNA in cellular compartments, 2) The distribution of mRNA along the skeleton of dendrites, 3) The degree of colocalization between two types of mRNA (not used in this paper) and 4) The degree of colocalization of mRNA at synapses in the dendrites.

Program Layout

The program is composed of 6 modules, two of which are executable, part1.py and part2.py. Both Part 1 and Part 2 draw from ui.py and segmentation.py. The following graph describes the purpose of each module and the module dependencies.




Figure 1. Schematic of the program layout, the function of each module and their dependencies. Only part1.py and part2.py are executable

Part 1

Terminology:

A **print** is a 2D binary image for a specified cellular compartment, where there is a white pixel where the cellular exists and a black pixel everywhere else

A **skeleton** is a 2D binary image of a single white line representing the midline of a print for a dendrite

MAP2 is an immunofluorescence using an antibody that detects the Microtubule associate protein 2 (Map2) which localizes in dendrites

DAPI (4',6-diamidino-2-phenylindole) is a blue-fluorescent DNA stain used to identify nuclei

Generating prints:

For dendrites, the program will look at each color used to annotate a dendrite in the annotation image and uses that color block as a mask for the 2D projection of the MAP2 image. In other words, the program extracts the map2 signal colored over by the annotation, and it does so separately for each dendrite. This MAP2 signal is subsequently binarized using the Otsu threshold generated from the original unmasked 2D MAP2 image. Thus, we obtain a binary image for each dendrite which is cleaned, smoothed, and then passed through an algorithm to account for inconsistent MAP2 signal.

Since synapses protrude out of dendrites, often the synapses are not contained by the Map2 stain, and thus also not contained by the dendrite prints. To account for this, when analyzing synapses, Part 1 of our program generates two prints for each dendrite, a regular one and another which is slightly dilated so that it contains all synaptic protrusions.

Because somas are easy to annotate accurately, when generating prints for somas, the program simply uses the masks generated from the annotations as a print. Thus, no binarizing is necessary, and only minimal cleaning is required to smooth the edges of the print.

For nuclei, the program will look at each color used to annotate a soma as a mask on the 2D projection of the DAPI channel. So, the program isolates the nuclei associated with each soma. The masked image is then binarized using a scaled Otsu threshold of the 2D DAPI image. Qualitatively, we found that scaling the threshold down by 0.65 worked best. DAPI stains tend to be stronger around the perimeter of the nucleus and spotty in the center, so once a binary image of each nucleus is obtained, it is cleaned, by filling in holes and smoothing edges.




Figure 2. Schematic of input and output of part1.py. Images of annotated dendrites and somas are used to generate prints and skeletons.

Algorithm for inconsistent signal

The MAP2 stain for dendrites varies in intensity even along the same dendrite. To obtain connected prints for each dendrite we merge fragments of the dendrite print with the following algorithm:

Input: 2D projection of a dendrite binarized by Otsu threshold

Algorithm:

- Save a copy of the input
- Obtain an over-dilated dendrite by dilating every pixel by a radius of 4 and then erode every pixel by a radius of 3 until all the white pixels in the image are part of 1 connected component
- Skeletonize the over-dilated dendrite
- Dilate the skeleton by a 5-pixel radius
- Add the dilated skeleton to the copy of the original input and obtain the final print




Figure 3. Schematic of the algorithm for inconsistent MAP2 signal

We choose to dilate by a radius of 4 pixels and erode by a radius of 3 pixel because it both smooths the edge of the image while also having the net effect of dilating each pixel in the dendrite by a radius of 1. Also note that because the dendrite is over-dilated by expanding each pixel equally in all directions, the midline (ie skeleton) of the original input should not be greatly affected by the algorithm.

Trimming Skeletons

Initial skeletons are obtained using the `skimage skeletonize` function. However, sometimes these original skeletons have some small branches. To trim these branches, use the `FilFinder` library to extract the part of the skeleton corresponding to the longest path. Additionally, the initial skeletons on occasion contain small circle or line segments disjoint from the main skeleton. In these cases, our `CleanBinaryImage` function is applied, which keeps only the largest component of the skeleton image.

Cleaning Binary Images

Even though every dendrite print is passed through the algorithm for inconsistent signal the final prints may still have white blobs which are disjoint from the main print. This can happen when the saved copy of the original dendrite print captures a signal that is noise and not from the dendrite. So even though the algorithm for inconsistent signal dilates the print into 1 connected component, when the skeleton is added back to the saved copy, the extraneous white blobs are still present.




Figure 4. Example noisy signal spots surrounding print after applying the algorithm for inconsistent signal

This issue is accounted for with a function for cleaning binary images, named `CleanBinaryImage`. This function looks at every connected component in the print and deletes all components except for the one with the largest area. However, it first checks that the second-largest component has an area smaller than 20% of the area of the largest component. If the print does not pass this requirement, it is evidence that there is an abnormally large extraneous signal or that something else has gone wrong in the generation of the print. In this case, the program will not delete any small blobs from the print, will report the problem to the user, and will not add this dendrite to any of the outline files.

Note the program also contains a similar function for refining soma and nucleus prints, named `refineSomaOrNucPrint` which calls the `Clean Binary Image` function. Because DAPI stains tend to be spotty at center of the nucleus, in addition to deleting small white blobs the program also fills in holes

Saving Skeletons and Prints

All prints and skeletons are saved as black and white gif images in a folder called `SkeletonsAndPrints`, so that they can be accessed by the second part of the program.

FISHQuant

FISHQuant is a program that searches smFISH images for Gaussian distributions of signals to detect mRNAs (Mueller et al., 2013). Furthermore, the program can distinguish between multiple mRNA spots if their signals overlap, making it particularly useful for counting the number of mRNAs in transcription sites. We use version one of FISHQuant, which is written in Matlab, to find the coordinates of all the mRNA in each dendrite and soma.

Part 2

This part of the pipeline can be used to extract different information and statistics from the mRNA/synapse coordinates and dendrite morphology to provide insights for different biological questions.

Currently functionalities include finding:

1. Density of mRNAs within each dendrite and soma under different experimental conditions
2. Distribution of mRNAs along the dendrites based on how far they travel from the soma
3. Colocalization of different mRNA species based on how far each mRNA is from its closest mRNA from the other species
4. Colocalization of mRNAs with synapses based on how many mRNAs are within a threshold distance from each synapse

Input Extraction

We first parse the information in the text files generated by performing spot detect in FISHQuant which records the coordinates of every mRNA and synapse. Thus, we first extract all the coordinates and which dendrite or soma they came from. It should be noted that the output of FISHQuant's spot detection function flips the x and y coordinates. We account for this in our own program.

```
FISH-QUANT
File-version      3D_v1
RESULTS OF SPOT DETECTION PERFORMED ON 17-Aug-2021
COMMENT Automated out-line definition (batch or quick-save)
IMG_Raw           21-08-07_DIV14_MG132_HSPABHSP110MAP2Tau_001_xy6_Cy5.tif
IMG_Filtered      21-08-07_DIV14_MG132_HSPABHSP110MAP2Tau_001_xy6_Cy5_filtered_batch.tif
IMG_DAPI          21-08-07_DIV14_MG132_HSPABHSP110MAP2Tau_001_xy6_DAPI.tif
IMG_TS_label      ..
FILE_settings     .._FQ_batch_settings_MATURE_2108017.txt
PARAMETERS
Pix-XY            300
107.5
CELL_START        Cell_1
X_POS            1146 1136 1135 1136 1142 1148 1152 1156 1160 1164 1167 1170 1177 1182 1186 1191 1196 1200 1205
1208 1211 1215 1219 1223 1226 1230 1233 1237 1241 1246 1251 1255 1259 1263 1267 1272 1279 1286 1293 1299
1304 1309 1315 1321 1327 1334 1340 1348 1357 1362 1368 1377 1385 1386 1385 1384 1383 1381 1380 1381 1382
1384 1387 1391 1394 1397 1401 1406 1412 1422 1431 1430 1426 1422 1420 1416 1413 1410 1408 1407 1407 1409
1410 1411 1412 1412 1408 1400 1392 1395 1391 1381 1371 1365 1359 1355 1351 1345 1336 1330 1325 1320
1313 1306 1299 1293 1289 1285 1281 1277 1272 1267 1262 1258 1255 1251 1248 1245 1241 1237 1233 1227 1223
1218 1216 1212 1208 1203 1199 1196 1193 1189 1186 1182 1177 1171 1163 1157 1153 1149 1499 1489 1479 1469 1459 1449
Y_POS            1632 1628 1618 1608 1598 1589 1579 1569 1559 1549 1539 1529 1519 1509 1499 1489 1479 1469 1459 1449
1439 1429 1419 1409 1399 1389 1379 1370 1360 1350 1340 1330 1320 1310 1300 1290 1280 1270 1260 1250 1240
1230 1220 1210 1200 1190 1180 1170 1160 1150 1140 1130 1120 1110 1100 1090 1080 1070 1060 1050 1040 1030
1020 1010 1000 990 980 970 960 950 946 954 964 974 984 994 1004 1014 1024 1034 1044 1054 1064
1074 1084 1094 1104 1114 1124 1134 1144 1154 1164 1171 1178 1186 1198 1208 1218 1226 1238 1248 1258
1268 1278 1288 1298 1308 1318 1328 1338 1348 1358 1368 1378 1388 1398 1408 1418 1428 1438 1447 1457 1467
1477 1487 1497 1507 1517 1527 1537 1547 1557 1567 1577 1587 1597 1607 1617
Z_POS
CELL_END
SPOTS_START
Pos_Y Pos_X Pos_Z AMP BGD RES SigmaX SigmaY SigmaZ Cent_Y Cent_X Cent_Z MuY MuX MuZ ITER_Y_det Y_det X_det Z_det Y_min
181608 152534 2080.29 1855.32 2167.74 1.27720e+07 145.371 145.171 469.374 212.892 217.42 935.874 209.894 186.569 1108.29 19 947 1420 10 945 949
1418 1422 7 13 4048 1568 384.017 0.78652 1 1 -1 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
101842 152888 5366.55 1144.3 1807.19 3.98548e+06 177.385 177.385 734.614 288.394 217.423 863.847 147.054 230.446 566.551 13 949 1423 20 947 951
1421 1425 17 23 3152 1083 211.098 0.398037 1 1 -1 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
102228 152791 4380.86 1583.02 1500.37 5.52946e+06 176.124 176.124 629.87 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
1415 1419 12 10 3124 1357 282.687 0.528092 1 1 -1 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
102307 152724 4466.42 1045.38 1932.88 4.18547e+06 176.124 176.124 629.87 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
1420 1424 13 19 2972 919 172.12 0.316669 1 1 -1 212.796 227.612 910.361 210.421 285.004 1000.86 13 952 1417 15 958 954
102645 152730 5324.5 1358.95 1150.66 5.64583e+06 221.167 221.167 868.026 211.809 227.703 885.021 197.023 294.936 824.498 17 956 1417 19 954 958
1415 1419 16 22 3038 1399 243.953 0.448829 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
102651 152757 5882.52 1383.38 1931.89 3.69842e+06 99.8517 99.8517 421.017 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
1420 1424 17 23 2973 1020 210.634 0.402246 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
102758 153186 5877.13 1181.59 1826.66 2.81195e+06 162.434 162.434 455.275 223.315 214.826 901.768 310.716 213.25 877.132 9 956 1426 18 954 958
1425 1429 15 21 3865 998 106.358 0.361263 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
102984 153732 3039.35 1490.67 2173.84 8.49488e+06 142.274 142.274 976.877 216.2 228.433 892.006 213.814 329.784 939.351 10 959 1430 14 957 961
1428 1432 11 17 3817 1353 369.086 0.680558 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
102984 153697 3562.77 1717.4 1477.84 5.96186e+06 266.892 266.892 1111.04 203.885 221.875 978.115 106.338 294.218 1762.77 12 960 1430 10 958 962
1428 1432 7 13 3262 962 100.933 0.55166 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
103283 153245 5912.66 1295.31 1878.27 1.31519e+07 199.973 199.973 583.651 209.4 209.633 926.082 190.392 165.103 1112.66 15 962 1427 20 960 964
1425 1429 17 23 3486 1346 354.786 0.65274 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
103452 153188 3813.92 1331.32 952.534 2.81709e+06 162.15 162.15 742.932 222.873 236.848 884.348 251.644 313.366 813.919 13 963 1413 14 961 965
1411 1415 11 17 2367 1068 223.562 0.411313 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
103710 152359 5755.58 1183.05 1296.57 2.49922e+06 154.334 154.334 677.837 210.296 220.527 863.1 187.068 246.673 655.579 14 966 1418 21 964 968
1416 1420 18 24 2573 1152 210.045 0.386443 1 1 -1 211.716 217.407 900.037 203.269 214.236 1082.52 12 956 1422 20 954 958
103850 153193 3816.45 1359.84 1269.04 4.57739e+06 149.311 149.311 596.772 216.332 222.509 903.197 219.989 250.038 916.452 11 967 1414 11 965 969
```

Figure 5. Example of output file from FISHQuant. The x and y coordinates of the mRNA or synapse signals are extracted for analysis.

Calculating Density

To calculate the area of a dendrite or soma, the program counts the number of pixels in the corresponding print generated from part 1. The area in pixels can then be converted to squared nanometers (nm^2) based on the conversion factor specific to each microscope. We then can take the ratio of the number of mRNA in and the area of each dendrite/soma to obtain compartment-specific mRNA density.

Number	Channel	Number of mRNA	Area (sq. nanometer)	Density
3	CY3	36	2808168.75	0.0000128197424033011
3	CY5	36	2808168.75	0.0000128197424033011
4	CY3	36	2657937.5	0.0000135443365391398
4	CY5	36	2657937.5	0.0000135443365391398
6	CY3	36	2877506.25	0.0000125108329477999
6	CY5	36	2877506.25	0.0000125108329477999

Figure 6. Example of compartment-specific mRNA count, area (in pixels) and mRNA density.

Finding the Endpoints of Dendritic Skeleton

To find the two endpoints of the skeleton, the program looks at the neighbors of each white pixel in the skeleton. If a skeleton point is an endpoint, only one of its neighboring pixels will be white, while all other points in the skeleton will have two neighboring white pixels. Since the skeleton can have a curved structure, the neighboring points of a skeleton point can be anywhere inside a 3 by 3 grid centered at the skeleton point of interest. As a result, the program counts the number of white pixels in a 3 by 3 grid centered at each point in the skeleton. If there are 3 white pixels in the grid, the point of interest is considered to not be an endpoint. If there are only 2 white pixels, the skeleton point centered is an endpoint.

If the program finds less than 2 or more than 2 endpoints for a particular skeleton it is an indicator that the outline for this dendrite was not generated properly. In this case, the program will not include the dendrite associated with the skeleton in the statistical analysis.

Determine the Soma Endpoint

Once two endpoints of the skeleton are found, the program aims to determine which of the two endpoints belongs to the soma end. Since dendritic annotation does not keep track of the soma end, the program needs to determine it using available information. With thousands of images from previous smFISH experiments, it has been observed that the soma end of each dendrite almost always has more mRNAs than the distal end. Biologically, this observation can be explained by the fact that the

transportation of mRNAs out of the soma is limited. Thus, we assume that the end of the dendrite with more mRNA is the end closer to the soma.

Distance of mRNA to Soma

To find how far the mRNAs are from the soma, the program first projects each mRNA to the closest pixel on the skeleton. Then it counts the number of pixels along the skeleton from the soma endpoint to the projected mRNA. As we count pixels along the skeleton, we keep track if an adjacent or diagonal pixel step was taken, so that we can calculate the precise distance of each step ($\sqrt{2}$ for a diagonal step and 1, for an adjacent step). Note that the distance along the skeleton from an mRNA to an endpoint is more accurate than the Euclidean distance because it considers the curved morphology of the dendrite. Since the skeleton is the center axis of the dendrite, the skeleton and dendrite share the same curvature. The calculated distances are then binned into groups of 25 micrometers before being written to an excel file.




Figure 7. A) Diagram of mapping mRNA spot coordinates to the skeleton B) Measuring distance of mRNA to soma by counting pixels along the skeleton

Dendrite Num	Channel	0-25 (um)	25-50 (um)	50-75 (um)	75-100 (um)	100-125 (um)	125-150 (um)	>= 150 (um)
1	Cy3	89	35	33	32	35	22	0
1	Cy5	47	9	8	13	10	6	0
2	Cy3	52	36	42	2	0	0	0
2	Cy5	45	16	12	0	0	0	0
3	Cy3	67	40	26	48	2	0	0
3	Cy5	61	25	13	25	0	0	0
4	Cy3	110	0	0	0	0	0	0
4	Cy5	68	0	0	0	0	0	0
5	Cy3	90	64	68	71	55	0	0
5	Cy5	51	40	26	35	32	0	0
6	Cy3	74	43	34	30	31	16	0
6	Cy5	38	14	13	4	9	3	0
7	Cy3	59	36	35	0	0	0	0
7	Cy5	52	20	10	0	0	0	0
8	Cy3	181	0	0	0	0	0	0
8	Cy5	96	0	0	0	0	0	0
9	Cy3	101	69	49	30	4	0	0
9	Cy5	73	26	12	5	1	0	0

Figure 8. Example of mRNA distribution based on distance to soma along the skeleton, each dendrite is assigned a unique number

mRNA Colocalization Statistics

Our motivation for calculating colocalization statistics for two kinds of mRNA is to see if the mRNAs are moving together down the dendrites in granules. The program will loop through every mRNA of species A and find the distance to the closest mRNA of species B. The same calculation can be done the other way around, finding the distance between each mRNA of species B to the closest mRNA of species A. Additionally, the program calculates self-colocalization, the distance between an mRNA of species A to another mRNA of species A (not including itself). To have a computational control, the same colocalization statistics were calculated with simulated mRNA coordinates. The simulation process is described below in the section "Simulation as Control." Nearest distances are then binned into user-defined groups before being written to an excel file.

Distance (nm)	Cy3 to closest Cy5	Cy5 to closest Cy3	Cy3 to closest Cy3	Cy5 to closest Cy5	Sim-Cy3 to closest Sim-Cy5	Sim-Cy5 to closest Sim-Cy3
0-75	487	489	94	109	0	0
75-150	1172	1190	387	359	0	0
150-225	1216	1198	710	550	0	0
225-300	1203	1113	988	671	0	0
300-375	792	655	1183	628	0	0
375-450	516	308	1191	643	0	0
450-525	448	187	1108	587	0	0
525-600	385	107	1087	518	0	0
600-inf	2863	277	2910	2168	9082	5524

Figure 9. Example of mRNA colocalization statistics: mRNAs are binned based on Euclidean distance to the closest mRNA of the other species of interest

Synapse Statistics

We seek to understand how often synapses are being served by mRNA. We assume that mRNA localizing near a synapse are servicing that synapse. We take advantage of the fact that PSD95 staining for a single synapse looks like a single mRNA spot and detect the position of synapses using FISHquant. When part2.py is run, the program will ask the user to define a colocalization threshold which will be used to determine whether a certain mRNA is localized at a synapse. Then the program loops through every synapse coordinate of each dendrite and counts the number of mRNA which are within a radius of the threshold distance. We record the number of synapses with X localized mRNA, for X equal to 0, 1, 2, 3, etc. As with the mRNA colocalization functionality, a computation control is required to interpret the statistical significance of the degree of synaptic localization. We generate fake mRNA coordinates, and again calculate the number of synapses being served by X fake mRNA.

Num of mRNA	Real Cy3	Sim Cy3	Num of mRNA	Real Cy5	Sim Cy5
0	1724	2164.61	0	2148	2394
1	525	109.73	1	358	140.48
2	150	78.5	2	100	60.24
3	32	46.14	3	16	23.12
4	5	23.15	4	6	7.87
5	1	10.21	5	0	2.4
6	1	3.84	6	1	0.64
7	0	1.39	7	0	0.2
8	0	0.32	8	0	0.04
9	0	0.1	9	0	0.01
10	0	0			
11	0	0			
12	0	0.01			

Figure 10. Example of synapse colocalization statistics: synapses are binned based on how many real or simulated mRNAs of interest are within its threshold radius

Simulation as Control

It is necessary to verify that the colocalization patterns (between mRNA and to synapses) observed are due to meaningful biological mechanisms (e.g., active transport of mRNA to synapse, co-transportation of different mRNA species) rather than restricted space and randomness. To do this, the program picks N fake mRNA sub-pixel coordinates uniformly random from each print, where N is the number of real mRNA in that dendrite. The program then uses these fake mRNA coordinates to calculate statistics in the same way it would with real mRNA coordinates. This simulation is repeated 100 times (50 times in the older version of the program), and averages of each statistic are taken. The simulation serves as computation control because if the simulated mRNAs show the same degree of colocalization as the real mRNAs, then colocalization patterns are likely due to having many mRNAs crammed into a restricted space. However, if the real mRNAs show strong colocalization, while simulated mRNAs do not, it suggests the existence of an underlying biological mechanism supporting mRNA localization.




Figure 11. Steps used to uniformly randomly sample sub-pixel coordinates from dendrite print to use as simulated mRNA spots

Writing to Excel Files

The results of all distribution and colocalization analysis are organized into tables and then saved as excel files. The tabular results can then be used to generate plots for visualization and examination.

Disclaimer: The program layout below is for the newer modularized version of our code. An older version of the code was used to collect the data presented in the paper. Both versions used the same fundamental algorithms, but the older version was organized so that each functionality of part2 was in its own executable script.

References

Mueller, F., Senecal, A., Tantale, K. *et al.* FISH-quant: automatic counting of transcripts in 3D FISH images. *Nat Methods* 10, 277–278 (2013). <https://doi.org/10.1038/nmeth.2406>

